

Tel: +1-262-636-9755
www.vistadental.com

VA
VISTA | APEX

 **Inter-Med, Inc.**
2200 South St
Racine, WI 53404 USA
Tel: +1-262-636-9755
www.vistadental.com

RX ONLY



Seamfree

Distributed by:
Inter-Med, Inc.
2200 South Street
Racine, WI 53404

VA
VISTA | APEX

Seamfree Wetting Resin and Lubricant

GENERAL INFORMATION

Seamfree is a wetting resin that has been designed to make the delivery of composite materials easier and more consistent by eliminating the tackiness that makes manipulation difficult. Not only does Seamfree eliminate handling issues but it doesn't build film thickness or add a visible layer between composite layers. Seamfree is compatible with all methacrylate materials and can be used in all direct and indirect restorations.

- Eliminates tackiness
- Adds no thickness or seams to the final restoration
- Compatible with all methacrylate materials
- For use in any restoration

PRINCIPAL USE OF Seam free :

- Lubricating restorative instruments and materials

Application Instructions:

Syringe:

1. Install applicator tip on the Seamfree syringe.
2. Dispense a small amount onto patient napkin to be sure material is flowing freely from the syringe.
3. Apply Seamfree to the desired surfaces to be lubricated or "de-tackified". The surface may be an instrument or uncured composite material.
4. Continue with the placement of the composite material. Activate all composite materials according to manufacturer's instructions. Seamfree will co-cure with the composite material, therefore, separate activation of Seamfree is not required.

Bottle:

1. Dispense one drop of Seamfree into a well.
2. Apply Seamfree to the desired surfaces to be lubricated or "de-tackified". An instrument may be dipped into the Seamfree or the Seamfree may be applied to an instrument or uncured composite material.
3. Continue with the placement of the composite material. Activate all composite materials according to manufacturer's instructions. Seamfree will co-cure with the composite material, therefore, separate activation of Seamfree is not required.

Note: Only a thin layer of Seamfree is necessary. Avoid excess build up to ensure acceptable results.

Inter-Med, Inc. recognizes its responsibility to replace products if proven to be defective. Inter-Med, Inc. does not accept liability for any damage or loss, either direct or consequential, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Revision 1: October 08, 2020

Safety Data Sheet

Section 1. Identification of Substance/Preparation and Company/Undertaking

Product Name: Seamfree
Manufacturer: Inter-Med, Inc.
2200 South Street
Racine, WI 53404

Emergency Telephone Number: 1-877-418-4782
Telephone Number for Information: 1-877-418-4782

Date Prepared: February 6, 2019

Section 2. Composition/Information on Ingredients

Ingredient	CAS No.	OSHA PEL	ACGIH TLV	Other
2-Methyl-2-Propenoic Acid (1 Methyl-2-ethylidene) Bis (4,1-Phenyleneoxy (2-Hydroxy-3,1-propanediyl)) Ester	1565-94-2	N/A	N/A	N/A
Triethylene Glycol Dimethacrylate	109-16-0	N/A	N/A	N/A

Section 3. Hazards Identification

Eye irritant
Possible skin irritant

Section 4. First-aid Measures

Inhalation: Provide fresh air and rest
Skin Contact: Flush with copious amounts of water
Eye Contact: Flush with running water for a minimum of 15 minutes. Seek medical attention
Ingestion: Drink 2 or more glasses of water or milk, seek medical attention. Do not induce vomiting

Section 5. Fire-fighting Measures

Fire Fighting Media: Dry chemical powder, foam or carbon dioxide
Fire and Explosion Risks: Upon combustion, irritating gases may form

Section 6. Accidental Release Measures

Precautions: Prevent skin/eye contact
Method for Clean-up: Wipe with damp cloth, discard appropriately

Section 7. Handling and Storage

Use according to directions for use. Store at room temperature in the original container provided by the manufacturer. Avoid unnecessary exposure to the product.

Section 8. Exposure Controls/Personal Protection

Eye Protection: Wear protective eyewear
Hand Protection: Wear protective gloves
Respiratory Protection: Avoid prolonged breathing of vapors

Skin Protection: Avoid contact, if contact does occur rinse with copious amounts of water

Section 9. Physical/Chemical Properties

Appearance and Odor: Tooth colored thin paste material
Boiling Point: > 100° Celsius
Vapor Pressure: (mm Hg) N/A
Vapor Density: (Air = 1) N/A
Specific Gravity: (Water = 1) 1
Melting Point: < 0° Celsius
Evaporation Rate: N/A
Solubility in Water: Not water soluble

Section 10. Stability and Reactivity

Stable. Avoid high temperatures

Section 11. Toxicological Information

Eye Contact: Irritates the eyes
Skin Contact: May irritate the skin
Sensitization: May cause allergic reaction upon contact with the skin
Other: N/A

Section 12. Ecological Information

Use as directed. Follow local law governing disposal of chemicals.

Section 13. Disposal Considerations

Follow local law and regulations for disposal of this material. If none exist follow state or regional laws and regulations

Section 14. Transport Information

Follow local law and regulations for disposal of this material. If none exist follow state or regional laws and regulations

Section 15. Regulatory Information

Classification: R36 Irritating to eyes., R38 Irritating to skin., R43 May cause sensitization by skin contact., Xi; Irritant,
Safety: S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. S45 In the event of an accident call immediately a physician. If possible show him/her this document.

Section 16. Other Information

Revision Date: January 28, 2020
Revision Number: 4

NOTE: This information is believed to be correct as of the date of this SDS. The manufacturer makes no representation as to the completeness or accuracy of this information and supplies it on the condition that the person receiving it will make their own determination as to its suitability or their purposes prior to use. In no event will the manufacturer be responsible for damages of any nature whatsoever resulting from the use of or reliance upon the information contained herein.