


exquisite

UNIVERSAL COMPOSITE

exquisite

UNIVERSAL COMPOSITE

REF 90496, 90497, 90498, 90499, 90500, 90501, 90502, 90503, 90504,
90505, 90506, 90507, 90508, 90509, 90510, 90511, 90580, 90581,
90582, 90583, 90584, 90585, 90586



Rx Only



Inter-Med, Inc.
2200 South St.
Racine, WI 53404, USA
Tel: +1-262-636-9755
www.vistaapex.com

NOTE: visit vistaapex.com for SDS

90386-I-AP-ENG-(5)

Description:

Exquisite Restoration® is a visible light-cure nanofilled composite for anterior and posterior restorations. Beside its low shrink and low abrasion Exquisite Restoration® has improved handling and features excellent mechanical properties. Exquisite Restoration® is highly radiopaque and is very wear resistant.

Indications for Use:

Exquisite Restoration® is suited for all Classes of direct restorations.

General Instructions¹:

Composite Fillings:

1. Isolate and prepare the tooth surface to be restored. Insert a proper matrix where indicated.
2. Use a cavity cleanser such as Clean & Boost® to make the cavity as clean as possible.
3. Apply an adhesive, such as Surpass® Self-Etch Adhesive system or Superb Universal Adhesive, according to the manufacturer's instructions.
4. In Class I's and Class II's insert a thin layer (< 0.5 mm) of flowable composite (such as Titan Flowable Composite) over the entire pulpal floor up to and including the cavosurface margins. Activate the flowable composite according to the manufacturer's instructions.
5. Insert a 2 mm increment of the desired shade of Exquisite Restoration® and properly adapt to the cavity.
6. Light-cure for 20 seconds.
7. Continue inserting increments and light curing until cavity is slightly over-filled.
8. Finish and polish as usual.

Storage:

Protect from sunlight and excessive temperatures (30°C/86°F). Store between 2°C/36°F and 30°C/86°F. Refrigerate materials that are not to be used for more than 7 days. See individual component labels for specific expiration dates. Do not use after the expiration dates. Close syringes immediately after use to avoid exposure to light.

Cautions:

Avoid contact with eyes. If the material comes into contact with the eyes, flush with copious amounts of water and seek medical attention. Un-polymerized resins may cause skin sensitization in susceptible persons, in case of contact with the skin, wash thoroughly with soap and water.

Notes:

1. Avoid the use of eugenol containing materials as eugenol can impair the hardening of resin materials and result in discoloration. Do not store resin based restoratives in the proximity of eugenol containing materials.
2. visit vistaapex.com for SDS information