

## **INSTRUCTIONS FOR USE: Please read completely before operating unit**

### **1. Warnings and Precautions**

- This system has been designed and manufactured to assure personal safety. Please read all safety and operating instructions carefully before installation and use.
- Handle system with care at all times.
- Follow normal operating procedures and asepsis techniques by using standard personal protection equipment such as: gloves, eyewear, facemasks and a protective gown while performing all dental procedures.

### **2. Unpacking the System**

- Gently unpack your Vista Syringe Warmer. Confirm that all components and accessories are included.

QTY	ITEM
1	Syringe Warmer Housing
1	6 or 12cc syringe dome
1	90° white power cord
5	Warmer base inserts (Reorder # 404317)

### **3. Operating your Vista Syringe Warmer**

- Take power cord and insert male end firmly into electrical input on syringe warmer housing.
- Place a single warmer base insert into the housing. Place warmer dome into housing.
- Plug warmer into standard wall outlet.
- Turn system on by switching the green ON/OFF switch on the front of the Syringe Warmer base to the "On" position. The switch will illuminate when the unit is on.
- Fill syringes with solution of choice.
- Place syringes into openings on syringe dome.
- Warmer will regulate temperatures between 110° and 130° F.

*Note:* Allow up to 20 minutes for system to reach optimal temperature when turning it on at the beginning of the day.

### **4. Indications**

- The Vista Syringe Warmer is indicated for the warm delivery of irrigating solutions.

### **5. System Care**

- It is recommended that the system be turned off when not in use.
- It is recommended that you wipe the system down at the end of each day.
- Do NOT attempt to autoclave the base or dome of the syringe warmer!
- Sanitize or dispose of syringes and needle tips between each patient.

### **Technical Support**

For questions regarding your Vista Syringe Warmer please call Vista Dental Products toll free at (877) 418-4782 Monday thru Friday 8:00 AM to 5:00 PM (Central Standard Time).

## **VISTA DENTAL PRODUCTS TERMS AND CONDITIONS OF WARRANTY:**

The operator assumes all risk and liability for damages arising out of the improper use of Vista's product. In the event of a defect in material or workmanship, Vista's liability is limited, at Vista's option, to replacement of the defective product, or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Vista.

The Syringe Warmer is warranted to be free from defects under normal usage conditions for 6 months from its date of delivery. There is no warranty, expressed or implied, of merchantability or fitness. The manufacturer's sole obligation under this warranty is to opt to either repair or replace the defective part of the product. If service must be performed to correct a defect, then the manufacturer will provide the service at its factory according to the mutual agreement made in advance. The manufacturer and its distributors will not accept the return of product unless the return is authorized and shipped in accordance with the distributor's instructions. Contact the local representative of the distributor from which the product was purchased for shipping instructions, a return authorization number, and an ARS shipping label.

There is no warranty, remedy or condition, expressed or implied, except as provided herein. The warranty and remedies contained herein are made by the manufacturer to the first buyer for dental use and are in lieu of all other agreements (expressed or implied), liabilities or remedies for breach of warranty. Vista Dental Products shall not be liable for consequential or incidental damages. No person or distributor is authorized to modify the terms of this warranty.

This warranty is void if any defect is caused by conditions beyond the manufacturer's control, including acts of God, damage resulting from mishandling, neglect, misuse, improper maintenance, accident or alteration/repair by anyone other than the manufacturer. The buyer assumes all liability for any damage caused by improper use of the product. The manufacturer assumes no liability for the user's failure to follow the instructions contained in this manual.

### **RETURN POLICY**

Vista Dental Products will accept for return previously purchased merchandise which is suitable for resale or that which was shipped in error by Vista Dental Products. Merchandise suitable for resale requires current labeling and unopened non-soiled packaging.

All returns must have prior approval and must be shipped "prepaid" along with a return authorization form and a copy of the original invoice. Any products returned that are discontinued, dated, damaged, or opened could be denied credit or assessed a higher return fee.

Merchandise returned for credit must be received by Vista Dental Products within 60 days of the original invoice date. Returns made within 30 days will be subject to a 15% restocking fee. Any returns made 31-60 days after the original invoice date will be subject to a 25% restocking fee.

Equipment cannot be returned without written authorization from Vista Dental Products. Any equipment returned within 30 days from the date of the original shipment from Vista may not be assessed a restocking fee as long as the merchandise has current labeling and unopened non soiled packaging. Unopened equipment returned within 31 to 60 days from the date of the original shipment from Vista requires a restocking fee of 25% of the purchase price, including shipping and handling charges. Any equipment returned after 60 days from date of original shipment from Vista will not be restockable for credit. Installation, if required, must be initiated with an outside party, and is the sole responsibility of the customer.

- Special orders are not suitable for resale and therefore not returnable for credit.
- Claims for lost or damaged shipments should be filed immediately with the carrier.