



VISTA | APEX

# Therma-Flo™ Composite Warming Kit



404350, 404355, 404310, 404320

**Rx ONLY**

## INSTRUCTIONS FOR USE:

**Please read completely before operating unit**

### 1. Warnings and Precautions

- This system has been designed and manufactured to assure personal safety. Please read all safety and operating instructions carefully before installation and use.
- Handle system with care at all times.
- The Therma-Flo™ Composite Warming kit operates at a high temperature and may be hot to the touch. Caution should be taken to avoid touching hot surfaces.

### 2. Indications

The Therma-Flo™ Composite Warming unit is designed to heat composite material used for dental restorations, to 155°F (68°C). Scientific literature and in-house testing shows heating of composites provides better flowability and stronger, longer lasting restoration in less time. There is a significant increase in flowability of highly filled composites through the use of the Therma-Flo™ kit. The warming unit accepts most composite syringes and the Step Down™ Tips will self tap on to most standard highly filled composite capsules.

### 3. Operating your Therma-Flo™ Composite Warming Kit

- Gently unpack the contents of your Therma-Flo™ Warming Kit.
- Confirm that all components and accessories are included:

#### Therma-Flo™ Complete Kit:

- (1) Composite Warmer Unit, (1) 90° white power cord, (2) Snap Fit Composite Syringes, (40) 4X Therma-Flo™ Step Down™ Tips

#### Therma-Flo™ Base:

- (1) Composite Warmer Unit, (1) 90° white power cord
- Connect power cord into the power connector making sure that it is fully seated.
- Turn system on by switching the gree ON/OFF switch on the front of the Composite Warmer base to the 'On' position. The switch will illuminate when the unit is on.
- Heat composite warming base for a minimum of 30 minutes prior to inserting loaded composite syringes. The unit should be turned on at the beginning of the day and off at the end of the day.

#### Warming Composites:

1. Insert capsule into Snap Fit composite syringe.
2. Insert Snap Fit composite syringe into an open slot in the composite warmer dome. Up to 3 capsule syringes may be inserted into warming unit for multi layer restorations. Warming unit is designed to hold additional capsules in the side warming slots of dome.
3. Heat capsules in composite syringe for 5-10 minutes prior to completing the restoration.

#### Applying Optional Step Down™ Tips:


1. Prior to heating, turn tip clockwise onto capsule until adequately secured and tight.
2. Follow steps above for warming composites.
3. Discard tip after use.

### 4. System Care

- It is recommended that the unit is turned off and system is wiped down at the end of each day.
- Do NOT attempt to auclave the base of dome of warming unit.

### Technical Support

For questions regarding your Therma-Flo™ Composite Warming unit, please call Vista Dental Products toll free at (877) 418-4782 Monday thru Friday 8:00 AM to 5:00 PM (Central Standard Time).



## **VISTA DENTAL PRODUCTS TERMS AND CONDITIONS OF WARRANTY:**

The operator assumes all risk and liability for damages arising out of the improper use of Vista's product. In the event of a defect in material or workmanship, Vista Apex's liability is limited, at Vista Apex's option, to replacement of the defective product, or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Vista Apex.

The Therma-Flo™ Composite Warming Kit is warranted to be free from defects under normal usage conditions for 6 months from its date of delivery. There is no warranty, expressed or implied, of merchantability or fitness. The manufacturer's sole obligation under this warranty is to opt to either repair or replace the defective part of the product. If service must be performed to correct a defect, then the manufacturer will provide the service at its factory according to the mutual agreement made in advance. The manufacturer and its distributors will not accept the return of product unless the return is authorized and shipped in accordance with the distributor's instructions. Contact the local representative of the distributor from which the product was purchased for shipping instructions, a return authorization number, and an ARS shipping label.

There is no warranty, remedy or condition, expressed or implied, except as provided herein. The warranty and remedies contained herein are made by the manufacturer to the first buyer for dental use and are in lieu of all other agreements (expressed or implied), liabilities or remedies for breach of warranty. Vista Apex shall not be liable for consequential or incidental damages. No person or distributor is authorized to modify the terms of this warranty.

This warranty is void if any defect is caused by conditions beyond the manufacturer's control, including acts of God, damage resulting from mishandling, neglect, misuse, improper maintenance, accident or alteration/repair by anyone other than the manufacturer. The buyer assumes all liability for any damage caused by improper use of the product. The manufacturer assumes no liability for the user's failure to follow the instructions contained in this manual.

### **RETURN POLICY**

Vista Apex will accept for return previously purchased merchandise which is suitable for resale or that which was shipped in error by Vista Apex. Merchandise suitable for resale requires current labeling and unopened non-soiled packaging.

All returns must have prior approval and must be shipped "prepaid" along with a return authorization form and a copy of the original invoice. Any products returned that are discontinued, dated, damaged, or opened could be denied credit or assessed a higher return fee.

Merchandise returned for credit must be received by Vista Apex within 60 days of the original invoice date. Returns made within 30 days will be subject to a 15% restocking fee. Any returns made 31-60 days after the original invoice date will be subject to a 25% restocking fee.

Equipment cannot be returned without written authorization from Vista Apex. Any equipment returned within 30 days from the date of the original shipment from Vista Apex may not be assessed a restocking fee as long as the merchandise has current labeling and unopened non-soiled packaging. Unopened equipment returned within 31 to 60 days from the date of the original shipment from Vista Apex requires a restocking fee of 25% of the purchase price, including shipping and handling charges. Any equipment returned after 60 days from date of original shipment from Vista Apex will not be restockable for credit. Installation, if required, must be initiated with an outside party, and is the sole responsibility of the customer.

- Special orders are not suitable for resale and therefore not returnable for credit.
- Claims for lost or damaged shipments should be filed immediately with the carrier.