



Rx ONLY

REF 317018 - 3% 16oz. Bottle

REF 317008 - 6% 16oz. Bottle

REF 281406 - 3% Gallon

REF 281403 - 6% Gallon

Description:

Sodium Hypochlorite is an aqueous solution.

It is used for the irrigation and debridement of root canals during and after instrumentation.

Procedure:

1. Irrigate with sodium hypochlorite as needed for debridement during instrumentation.
2. Remove the smear layer with appropriate irrigant, such as, EDTA solution or Citric Acid 12.5%.
3. Irrigate again with sodium hypochlorite for further debridement and cleansing.

Instructions for Use:


1. Select desired size syringe and irrigation tip.
2. Fill syringe with sodium hypochlorite.
3. Attach irrigation tip securely to syringe.
4. Irrigate as desired.
5. Remove and dispose of syringe and irrigating tip.

Precautions:

1. Read SDS for material handling and safety information.
2. Avoid allowing sodium hypochlorite to be ingested, use proper evacuation.
3. Avoid spilling on colored fabrics and surfaces.

Storage:

- Store in a cool, dry area away from heat and direct sunlight.
- If desired, transfer to syringe for short term storage, up to 2 weeks.



DANGER

Skin corrosion. Serious eye damage. Very toxic to aquatic life. Avoid release to the environment. Wear protective gloves/eye protection/face protection/protective clothing. Harmful if swallowed. Wash thoroughly after handling and avoid any chemical contact with eyes, skin and clothing. Keep container tightly closed to avoid inhalation or accidental ingestion. Use with adequate ventilation when necessary. Follow all government regulations for waste disposal. Discharge to the environment must be avoided. **FIRST AID: IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **IF ON SKIN:** Take off contaminated clothing and shoes immediately. Wash off with plenty of soap and water. Consult a physician. **IF INHALED:** Move person to fresh air. **IF SWALLOWED:** Do not induce vomiting. Rinse mouth with water. **Specific treatment:** Immediately call a Poison Control Center or doctor/physician. **FIRE: IN CASE OF FIRE:** Evacuate all personnel. Use dry powder, water spray, fog, mist, foam, sand or carbon dioxide to extinguish fire.