


Irrigation. Simplified.





Sodium Hypochlorite effectively degrades organic tissue and kills bacteria in the root canal system. Pharmacy controlled, dental indicated NaOCl to ensure consistency and safety.

Color coded with corresponding syringes* for safety and convenience

*syringes sold separately

Luer-Lock Cap to minimize waste and mess



EDTA 17% Solution removes the smear-layer that creates a better seal with bio-active sealers. Our EDTA is specifically engineered with a lower pH for optimal performance.



SCAN NOW to download PDF for easy implementation!

VistaApex
SOLUTIONS
Value. Added.

Sodium Hypochlorite 3% and 6% Concentration



Key benefits:

- For irrigation and debridement during and after instrumentation
- Provides tissue dissolution and solubilization of organic debris
- Quality controlled NaOCl for endodontic indications

Available in 3% Concentration 16oz (480mL) bottle (SKU 317018) or 6% Concentration 16oz (480mL) bottle (SKU 317008)

EDTA 17% Solution Disodium Ethylenediaminetetraacetic



Key benefits:

- Removes smear layer and dentin mud
- Opens dentin tubules and enhances cleansing
- Specifically engineered with lower pH for optimal performance^{1,2}

Available in 16oz (480mL) bottle (SKU 317011) or 4oz (120mL) bottle (SKU 317001)

Standard Irrigation Protocol

Standard protocol for tissue dissolution, smear layer removal and pristine cleanliness



STEP 1: Irrigate

During instrumentation and the initial irrigation, use 3% or 6% Sodium Hypochlorite to speed tissue dissolution.



STEP 3: Irrigate

Irrigate with VistaApex 3% or 6% Sodium Hypochlorite.



STEP 2: Clean

Remove the smear layer with VistaApex's EDTA 17% Solution.



STEP 4: Rinse

Rinse with sterile water or alcohol before obturation.



IRRIGATE

Use the VIP Tip 17mm (SKU 315729) or 27mm (SKU 315829) during any irrigation step.



ACTIVATE

Activation of the irrigants in Steps 1 and 2 of both protocols for 30 seconds can increase irrigation efficacy significantly (SKU 407702).



SEAL

Employ Bioglass (SKU 504753) to promote bio-regenerative healing and use our Material Saving Tip (MST) (SKU 312411) reducing waste by 78%.

VistaApex Solutions

2200 South Street
Racine, WI 53404

Phone: (877) 857-8213

Email: info@vista-dental.com

Website: <http://vistaapex.com>

References:

1. Nikiforuk G, Sreenby L. (1953). J Dent Research, 32(6): 859-867.
2. Serper A, Calt S. (2002). J Endod, 28(7): 501-502.