



THE POWER TO CONFIDENTLY CURE  
IS HERE

ONLY *ONE* IS FIT TO CURE.  
UNLEASH THE MULTISPECTRAL  
POWER THAT IS VALIANT®



# FEATURES


- Cordless
- Induction Charging Base
- One Battery Pack
- Warranty (3 Years)
- Cool Touch Technology
- Multispectral
- Light Weight
- Three Curing Modes
- Transillumination Light

# TECHNICAL SPECIFICATIONS

Range of Light Output (nm)	395 - 500
Weight	72g
Curing Modes	3 (Standard, Ramp, Boost)
Lens Size	10mm
Construction	Aluminum Body
Accessibility	10°
Battery Power	> 150, 20 sec Cure Cycles

# BROAD-SPECTRUM CURING

Wavelength Spectrums of Valiant® Modes



# VALIANT<sup>®</sup>

*THE POWER TO CONFIDENTLY CURE.*






# VALIANT® VS. COMPETITORS

Handpiece	Price	Multispectral	Weight	LEDs	Curing Mode	Transilluminator
<b>VALIANT®</b>	<b>\$995.00**</b>	✓	72g	4 + 1	3	YES
Valo® Brand	\$1581.00*	✓	170g	4	3	NO
Bluephase® Brand	\$1367.00*	✓	120g	4	1	NO
Smartlite® Focus® Brand	\$1365.00*	NO	87g	1	1	NO
Demi™ Ultra Brand	\$1650.00*	NO	184g	1	1	NO

\*\*Plus a \$200.00 Vista Merchandise Credit (excludes equipment).

\*Ultradent® and Valo® are registered trademarks of ULTRADENT Products, Inc. \*Bluephase® Style is a registered trademark of Ivoclar Vivadent Inc. \*Smartlite® Focus® is a registered trademark of DENTSPLY International. \*Demi™ Ultra and Kerr™ are both registered trademarks of Kerr Corporation. \*All saving percentages are based on Vista's quantity discount prices and competitors' 2018 prices.

## LED BEAM PROFILE



## ORDERING INFORMATION

### [600100] VALIANT® Curing Light Starter Kit

- (1) Cordless Handpiece
- (1) Battery Pack
- (1) Charging Base
- (100) Protective Sleeves
- (1) Instructions for Use
- (1) Universal Plug Set
- (1) Curing Light Shield

### [601001] VALIANT® Curing Light Shield

### [601004] VALIANT® Battery Pack

### [601008] VALIANT® Protective Barrier Sleeves



(877) 418-4782 | [www.vista-dental.com](http://www.vista-dental.com) | [info@vista-dental.com](mailto:info@vista-dental.com)

