

Tips ON TIPS



Endodontic | Periodontal | Tissue Management | Universal



UNIVERSAL TIPS

Appli-Vac™ – 3/4" Bendable Tip

- Blunt end, dead-soft needle resists crimping

Designed for: Sealants, flowable composites, etchants, resins, primers



100 pk

- 16ga - 315116
- 18ga - 315118
- 19ga - 315119
- 20ga - 315120
- 22ga - 315122
- 25ga - 315125

Secure-Lock™ Appli-Vac™ 3/4" Bendable Tip

- Increased safety and ease of use

Designed for: Sealants, flowable composites, etchants, resins, primers



20 pk

- 19ga - 314019
- 20ga - 314020
- 22ga - 314022
- 25ga - 314025

80 pk

- 19ga - 314119
- 20ga - 314120
- 22ga - 314122
- 25ga - 314125

Pre-Bent Tip

- Blunt end pre-bent needle
- Opaque photo sensitive hub (22/25ga)

Designed for: Etchants, resins and flowable composite



100 pk

- 16ga - 315316
- 22ga - 315322
- 20ga - 315320
- 25ga - 315325

1/2" Bendable Applicator Tip

- Blunt end, dead-soft needle resists crimping

Designed for: Endo Irrigating Solutions



100 pk

- 22ga - 315222
- 25ga - 315225

Sol-u-Flo™ Applicator Tip

- 19ga tip with adjustable brush fiber length
- Curved, cannula for optimal entry

Designed for: Aqueous solutions such as Caries Indicator, Silane & Chlorhexidine



100 pk - 312200, 500 pk - 312205

Exacta-Flo™ Brush Tip

- Adjustable, soft brush fibers
- Minimizes bubbles

Designed for: Aqueous solutions such as Caries Indicator, Silane & Chlorhexidine



20 pk - 321620, 100 pk - 321610, 500 pk - 321650

Volu-Flo™ Tip

- Delivers large volumes of material (18ga)
- Photo sensitive tip maintains the flow of light-cured material

Designed for: Viscous materials, pastes & putty's



100 pk - 321100

Blue-Flo™ Tip

- 23ga soft plastic tip
- Autoclavable

Designed for: Irrigation, aeration and impression materials



20 pk - 316202, 100 pk - 316222

Spira-Flo™ Brush Tip

- Nylon brush fibers for controlled delivery
- Internal spiral chamber minimizes bubbles

Designed for: Etchants & resins



100 pk - 321400, 500 pk - 321450

TISSUE MANAGEMENT TIPS

Stat-Flo™ Tip

- 19ga tip
- Padded-end brush tip

Designed for: Hemostatic agents




20 pk - 312102, 100 pk - 312100, 500 pk - 312105

ENDODONTIC TIPS

NiTi SuperFlex™ Tip

- **AUTOCLAVABLE.** Slotted and side-vented, 30ga for safe irrigation.
- 2.5x more flexible

Designed for: Endo Irrigating Solutions



MST Visco-Tip™

- Bendable, flexible 25ga tip greatly reduces waste
- Improves flow and placement of viscous materials

Designed for: Calcium Hydroxide and Bioceramic materials



Elasti-Vac™ Tip

- Allows for quick and easy removal of moisture from endodontic canal systems
- Narrow, flexible tip


Designed for: Endodontic evacuation



Appli-Vac™ 1" Bendable Tip

- Slotted and side-vented
- Individually capped in 100pk

Designed for: Endo Irrigating Solutions



- | | |
|---------------|---------------|
| 23ga - 315023 | 23ga - 315523 |
| 27ga - 315027 | 27ga - 315527 |
| 30ga - 315030 | 30ga - 315530 |

Secure-Lock™ Appli-Vac™ 1" Bendable Tip

- Slotted and side-vented
- Individually capped in 80pk

Designed for: Endo Irrigating Solutions



- | | |
|---------------|---------------|
| 23ga - 314023 | 23ga - 314523 |
| 27ga - 314027 | 27ga - 314527 |
| 30ga - 314030 | 30ga - 314530 |

Flexi-Glide™ Utility Tip

- 30ga flexible tip for easy access in curved canals
- Tip can be cut to desired length

Designed for: Endo Irrigating Solutions



Vista-Probe™ 1" Bendable Tip

- Closed-end for side-port delivery
- Individually capped in 100pk

Designed for: Perio (23ga) and Endo Irrigating Solutions



- | | |
|---------------|---------------|
| 23ga - 315123 | 23ga - 315623 |
| 27ga - 315127 | 27ga - 315627 |
| 30ga - 315130 | 30ga - 315630 |

Secure Lock™ Vista-Probe™ 1" Bendable Tip

- Closed-end for side-port delivery
- Individually capped in 80pk

Designed for: Perio (23ga) and Endo Irrigating Solutions




- | | |
|---------------|---------------|
| 23ga - 314123 | 23ga - 314623 |
| 27ga - 314127 | 27ga - 314627 |
| 30ga - 314130 | 30ga - 314630 |

PERIODONTAL TIP

Pocket-Probe™

- Multiple perforations for thorough flushing of perio pockets
- 21ga, closed-end, flexible irrigation tip

Designed for: Oral Rinse solutions such as Chlorhexidine





MATERIAL SAVING TECHNOLOGY

Conventional tips trap a significant amount of unusable material inside of the hub.



Conventional tip



MST Visco-Tip™

- MST minimizes the internal volume of the tip greatly reducing waste, **up to 87%!**
- Materials go further, and you can use **more of what you paid for!**

CASE Study


How much savings will MST provide with popular materials?

Testing Method:

Comparative waste testing of Vista's MST tips vs. the standard tips provided by the material manufacturer. Tip mass was measured prior to, and following the extrusion of material. The difference in mass between these two measurements represents the amount of wasted material remaining in the tip.

PRODUCT	COST PER SYRINGE	MASS WASTED(mg) Standard vs MST	% SAVINGS
Brasseler®: BC Sealer®*	\$180	144/19 mg	87%
Bisco®: TheraCal®*	\$40	93/17 mg	82%
Pulpdent®: Embrace™*	\$22	57/15 mg	74%

ENDODONTIC TIPS




MST Visco-Tip™

- Bendable, flexible 25ga tip
- Improves flow and placement of viscous materials

Designed for: Bioceramic sealers, and Calcium Hydroxide

20 pk - 312302
50 pk - 312300

UNIVERSAL TIPS



MST Tip

- 22ga Pre-bent tip
- Opaque photo sensitive hub

Designed for: Costly flowable composites, pulp capping materials, etchants, sealants

100 pk - 315310

*Brasseler® is a registered trademark of Peter Brasseler Holdings, LLC. Bisco® and TheraCal® are registered trademarks of Bisco, Inc. Pulpdent® is a registered trademark of Pulpdent Corporation. Embrace™ is a trademark of Pulpdent Corporation.