

# BE SMOOTH BE TACK-FREE



**Seamfree**  
WETTING RESIN & LUBRICANT

- Eliminates tackiness of composite materials
- Does **NOT** add layers or thickness to composite
- Absorbs into composite after use, does **NOT** need to be dried or light activated
- Does **NOT** change physical properties of composite material

**VISTA | APEX**  
vistaapex.com | 877.418.4782

# 00000

**SeamFree** is a wetting resin that has been designed to make the handling of composite material easier and more consistent by eliminating the tackiness that makes manipulation difficult. Not only does SeamFree eliminate handling issues, but it doesn't build film thickness or add a visible layer between composite layers.

	<b>SeamFree</b>	<b>Ultradent®* Composite Wetting Resin</b>	<b>Bisco®* Modeling Resin</b>
<b>Filler Content</b>	Unfilled	45% filled	30% filled
<b>Contains Solvent</b>	No	No	No
<b>Indications for Use</b>	Direct / Indirect	Direct / Indirect	Direct / Indirect
<b>Compatible with all Methylmethacrylate Materials</b>	Yes	Yes	Yes
<b>Adds Film Thickness or Seams to Final Restoration</b>	No	Yes	Yes
<b>Method of Application</b>	Syringe with Brush Tip, Bottle with micro applicator brush	Syringe with Brush Tip	Syringe with Tip

*"In conjunction with the composite heater this makes the spreading and contouring of composite amazingly easier. If you don't have a composite warmer, SeamFree is a must, as it also cuts the amount of time needed to polish."*

– Dr. Michael Miyasaki

*"What makes SeamFree a genius product is that it will co-polymerize instantly when covered by a material that will light activate, but never polymerize on its own. So you never have a thick layer of unpolymerized resin, you spend less time sculpting and less time polishing"*

– Dr. John Flucke

*"If you haven't used SeamFree yet, use it TOMORROW. This stuff is great and will take your polishing and shaping time down significantly. You drop a little bit of SeamFree on your composite and you can shape however you want without dragging"*

– Dr. John Kanca



**SeamFree Syringe Kit**



**SeamFree Bottle Kit**

## Kit Options

90025-2

### SeamFree 1 Syringe Kit

- 1 x 1.2mL Prefilled Syringe
- 10 x Exacta-Flo™ Tips

90026-2

### SeamFree 5 Syringe Kit

- 5 x 1.2mL Prefilled Syringe
- 50 x Exacta-Flo™ Tips

90027-2

### 3mL Bottle Kit

- 1 x 3mL bottle
- 1 x Single Cavity Well
- 25 x Applicators

\*Ultradent® is a registered trademark of Ultradent Products, Inc., Bisco® is a registered trademark of Bisco, Inc.