

Equipment Eligible for **Section 179 Tax Deduction** - AND **15% OFF** through the End of the Year!

FEATURED PRODUCT



Pinkwave™ QuadWave™ Curing Light

PinkWave™ is the next generation of curing lights, equipped with patented QuadWave™ Technology, offering decreased composite shrinkage and increased depth of cure. With the largest curing light area on the market, PinkWave™ ensures maximum coverage. QuadWave™ Technology leverages four wavelengths into one collimated pink beam: UV, Blue, Red, and Near-infrared.

- 12.1mm curing area
- Decreased shrinkage
- Up to 8.5mm depth of cure

**15%
OFF**

PRODUCT	SKU
Pinkwave™ QuadWave™ Curing Light Kit	90734
<i>Includes: (1) Pinkwave™ Curing Light, (1) charging base, (5) autoclavable light shields, (100) disposable barriers</i>	



EndoUltra® Cordless Ultrasonic Activator

EndoUltra's patented technology produces ultrasonic tip frequency of up to 40,000 Hz - 200 times more powerful than sonic devices!

- Better cleansing in less time
- Better penetration of irrigants throughout the canal
- Better acoustic streaming with cavitation

**15%
OFF**

PRODUCT	SKU
EndoUltra® Kit	407702
<i>Includes: (1) EndoUltra® cordless handpiece, (3) Titanium activator tips - 20/02, (3) Autoclavable handpiece sleeves, (100) Protective barriers, (1) Tip wrench, (1) Charger</i>	



Phasor™ Composite Warming Gun

The Phasor™ utilizes near-infrared technology to rapidly warm composite compules!

- Warms to 150°F in under 1 minute
- Better adaptation
- Reduced voids and microleakage
- Improved depth of cure
- Compatible with most compules

**15%
OFF**

PRODUCT	SKU
Phasor™ Composite Warming Gun	408800
<i>Includes: (1) Phasor™ handpiece, (1) charging stand, (1) 2 pack of nosecones, (1) box of barrier sleeves, 100/box, (1) maintenance kit, (1) medical USB charger, (1) USB to Micro B cable</i>	

Equipment Eligible for **Section 179 Tax Deduction** - AND **15% OFF** through the End of the Year!



Therma-Flo™ Composite Warming Kit

This composite warmer was uniquely engineered to utilize heat for optimal performance of any preferred composite material.

15% OFF

PRODUCT	SKU
Therma-Flo™ Composite Warming Kit	404350

Includes: (1) Therma-Flo™ Composite Warming Unit (110V), (2) Composite Syringes, (40) 4X Step Down Tips

15% OFF



Syringe Warmer

The Vista Apex Syringe Warmer is perfect for heating solutions (up to 130°F) - ideal for irrigating solutions and anesthetics!

PRODUCT	SKU
6cc Syringe Warmer	404301
12cc Syringe Warmer	404302

Includes: (1) Syringe Warmer unit, (5) warmer bowl inserts

15% OFF

Fill Station™

This patented syringe pump fill station offers a simple, convenient way to fill a dental syringe with liquid solutions!



PRODUCT	SKU
Fill Station™	408400

Includes: (1) Fill Station™ Unit

15% OFF

Valiant® Multispectral Curing Light

- 395-500nm Wavelength
- 4 LEDs and 3 Curing Modes

PRODUCT	SKU
Valiant® Multispectral Curing Light	600100

Includes: (1) Valiant® cordless handpiece, (1) battery pack, (1) charging base, (100) protective sleeves

