

# COVID Essential Canal Clean™

## Two-Step Procedure Kit

"By increasing the usage of kits instead of buying individual components when possible, you eliminate the amount of handling of each individual product, reducing the risk of cross contamination. Look for kits like CanalClean™ for endodontic irrigation."

- Dr. Michael Miyasaki



## PEEL AND USE ...

- Roughly \$8 per procedure

## Steps toward simplified irrigation

- Single use procedural irrigation kits are a great way to increase confidence and efficiency in irrigation.
- CanalClean™ provides everything needed for safe, reliable irrigation in a peel and use pouch. The two step procedure featuring Chlor-XTRA™ and SmearOFF™ is simple to use, always fresh, and proven to be effective.
- Additionally, CanalClean™ eliminates prep-time, all of the solutions are pre-packaged in syringes, including needle tips, and even a micro-suction tip for drying the canals.
- A great solution for anyone looking to expand their endodontic offering.

**SKU QTY**

### Canal Clean Kit

[504700] 4 Kits/box

Each kit includes:

- (1) 12cc Chlor-XTRA™ (6% NaOCl), (1) 6cc SmearOFF™ (2-in-1),
- (2) 30ga side-vented tips, and (1) Micro Aspirator with tip

### Chlor-XTRA™ (patented)

- 3x thinner than standard NaOCl
- 2x greater digestion capability

### SmearOFF™ (patented)

- Combination EDTA/CHX
- Will not form a precipitate when mixed with NaOCl