



**THE POWER TO CONFIDENTLY CURE
IS HERE**

**UNLEASH THE MULTISPECTRAL
POWER THAT IS VALIANT**



VA VISTA | APEX

FEATURES


- Cordless
- Induction Charging Base
- One Battery Pack
- Light Weight
- Three Curing Modes
- Transillumination Light
- Cool Touch Technology
- Multispectral

TECHNICAL SPECIFICATIONS


Range of Light Output (nm)	395 - 500
Weight	72g
Curing Modes	3 (Standard, Ramp, Boost)
Lens Size	10mm
Construction	Aluminum Body
Accessibility	10°
Battery Power	> 150, 20 sec Cure Cycles

BROAD-SPECTRUM CURING

Wavelength Spectrums of Valiant® Modes




VALIANT[®]



***THE POWER
TO CONFIDENTLY CURE***

LED BEAM PROFILE



ORDERING INFORMATION

[600100-2] VALIANT® CURING LIGHT COMPLETE KIT

- (1) Cordless Handpiece
- (1) Battery Pack
- (1) Charging Base
- (100) Protective Sleeves
- (1) Instructions for use
- (1) Universal Plug Set
- (1) Curing Light Shield

[601001-2] VALIANT® CURING LIGHT SHIELD

[601004-2] VALIANT® BATTERY PACK

[601008-2] VALIANT® PROTECTIVE BARRIER SLEEVES