

# Tips ON TIPS



Endodontic | Periodontal | Tissue Management | Universal



## UNIVERSAL TIPS

### Appli-Vac™ – ¾" Bendable Tip

- Blunt end, dead-soft needle resists crimping

*Designed for: Sealants, flowable composites, etchants, resins, primers*



**100 pk**

- 16ga - 315116
- 18ga - 315118
- 19ga - 315119
- 20ga - 315120
- 22ga - 315122
- 25ga - 315125

### ½" Bendable Applicator Tip

- Blunt end, dead-soft needle resists crimping

*Designed for: Etchants, resins, and flowable composites*



**100 pk**

- 22ga - 315222
- 25ga - 315225

### Micro-Flo™ Tips

- For microscope enhanced dentistry

*Designed for: Etchants and composites in micro-restorative procedures*



**50 pk**

- 26ga - 406626    27.5ga - 406627    30ga - 406630

### Pre-Bent Tip

- Blunt end pre-bent needle

- Opaque photo sensitive hub (22/25ga)

*Designed for: Etchants, resins and flowable composite*



**100 pk**

- 16ga - 315316
- 20ga - 315320
- 22ga - 315322
- 25ga - 315325

### Volu-Flo™ Tip

- Delivers large volumes of material (18ga)
- Photo sensitive tip maintains the flow of light-cured material

*Designed for: Viscous materials, pastes & putty's*



**100 pk - 321100, 500 pk - 321650**

### Blue-Flo™ Tip

- 23ga soft plastic tip
- Autoclavable

*Designed for: Irrigation, aeration and impression materials*



**100 pk - 321100**

### Spira-Flo™ Brush Tip

- Nylon brush fibers for controlled delivery
- Internal spiral chamber minimizes bubbles

*Designed for: Etchants & resins*



**100 pk - 321400, 500 pk - 321450**

## TISSUE MANAGEMENT TIPS

### Stat-Flo™ Tip

- 19ga tip
- Padded-end brush tip

*Designed for: Hemostatic agents*



**20 pk - 312102, 100 pk - 312100, 500 pk - 312105**

## ENDODONTIC TIPS

### NiT SuperFlex™ Tip

- AUTOCLAVABLE. Slotted and side-vented, 30ga for safe irrigation.
- 2.5x more flexible

*Designed for: Endo Irrigating Solutions*



### MST Visco-Tip™

- Bendable, flexible 25ga tip greatly reduces waste
  - Improves flow and placement of viscous materials
- Designed for: Viscous materials such as bioceramic sealers and calcium hydroxide*



### Elasti-Vac™ Tip

- Allows for quick and easy removal of moisture from endodontic canal systems
- Narrow, flexible tip

*Designed for:*  
*Endodontic evacuation*



### Appli-Vac™ 1" Bendable Tip

- Slotted and side-vented
- Individually capped in 100pk

*Designed for: Endo Irrigating Solutions*



### Secure-Lock™ Appli-Vac™ 1" Bendable Tip

- Slotted and side-vented
- Individually capped in 80pk

*Designed for: Endo Irrigating Solutions*



### Flexi-Glide™ Utility Tip

- 30ga flexible tip for easy access in curved canals
- 17mm tip can be cut to desired length

*Designed for: Endo Irrigating Solutions*



### Vista-Probe™ 1" Bendable Tip

- Closed-end for side-port delivery
- Individually capped in 100pk

*Designed for: Perio (23ga) and Endo Irrigating Solutions*



### Secure Lock™ Vista-Probe™ 1" Bendable Tip

- Closed-end for side-port delivery
- Individually capped in 80pk

*Designed for: Perio (23ga) and Endo Irrigating Solutions*



## PERIODONTAL TIP

### Pocket-Probe™

- Multiple perforations for thorough flushing of perio pockets
- 21ga, closed-end, flexible irrigation tip

*Designed for: Oral Rinse solutions such as Chlorhexidine*





# MATERIAL SAVING TECHNOLOGY

Conventional tips trap a significant amount of unusable material inside of the hub.



Conventional tip



MST Visco-Tip™

- MST minimizes the internal volume of the tip greatly reducing waste, **up to 87%!**
- Materials go further, and you can use ***more of what you paid for!***

## CASE Study

*How much savings will MST provide with popular materials?*

### Testing Method:

Comparative waste testing of Vista's MST tips vs. the standard tips provided by the material manufacturer. Tip mass was measured prior to, and following the extrusion of material. The difference in mass between these two measurements represents the amount of wasted material remaining in the tip.

PRODUCT	COST PER SYRINGE	MASS WASTED(mg) Standard vs MST	% SAVINGS
Brasseler®: BC Sealer®*	\$180	144/19 mg	87%
Bisco®: TheraCal®*	\$40	93/17 mg	82%
Pulpdent®: Embrace™*	\$22	57/15 mg	74%

### ENDODONTIC TIPS



**MST Visco-Tip™**

- Bendable, flexible 25ga tip
- Improves flow and placement of viscous materials

*Designed for: Bioceramic sealers, and Calcium Hydroxide*

20 pk - 312302    50 pk - 312300

### UNIVERSAL TIPS



**MST Tip**

- 22ga Pre-bent tip
- Opaque photo sensitive hub

*Designed for: Costly flowable composites, pulp capping materials, etchants, sealants*

100 pk - 315310

\*Brasseler® is a registered trademark of Peter Brasseler Holdings, LLC. Bisco® and TheraCal® are registered trademarks of Bisco, Inc. Pulpdent® is a registered trademark of Pulpdent Corporation. Embrace™ is a trademark of Pulpdent Corporation.