



VISTA APEX
ENDO RESTO LET'S GO!



pinkwave™
led curing light

25% OFF

Complete Kit
Originally \$1399⁰⁰

\$1049²⁵

PinkWave™ is the next generation of curing lights, equipped with patented QuadWave™ Technology, offering decreased composite shrinkage and increased depth of cure. With the largest curing light area on the market, PinkWave™ ensures maximum coverage. QuadWave™ Technology leverages four wavelengths into one collimated pink beam: UV, Blue, Red, and Near-infrared.

PinkWave™ is the official curing light of SmileFast USA

PRODUCT	SKU	25% OFF PRICE
Complete Kit <i>Includes: (1) Curing Light, (1) Charging Base, (5) Autoclavable Light Shields, (1) Pack of 5 Anti-Roll Brackets, (100) Disposable Barriers and (3) Complementary Pairs of Safety Glasses</i>	90734	\$1049²⁵
PinkWave™ Paddle <i>Includes: (5) Autoclavable Light Shields</i>	90786	\$26²⁴
Disposable Barriers <i>Includes: (5) Autoclavable Light Shields</i>	90752	\$19⁸⁹
See website for more product options		

Quick-Stat™ FREE - Clear Hemostatic Gel



25% OFF

Standard Kit
Originally \$23⁴⁸
\$17⁶¹

Quick-Stat™ FREE is the first clear tissue management agent with shortened hemostasis time. Not only that, the proprietary surfactant in Quick-Stat™ FREE allows better penetration to more effectively stop sulcular fluid build up. This is the Astringent Hybrid Technology that makes Quick-Stat™ FREE 20% faster than other clear hemostatic agents. For controlling tissue management, especially in the aesthetic zone, Quick-Stat™ FREE is the clear choice.

PRODUCT	SKU	25% OFF PRICE
Standard Kit <i>Includes: (4) 1.2mL Prefilled Syringes and (8) Stat-Flo™ Tips</i>	504600	\$17⁶¹
Value Pack <i>Includes: (40) 1.2mL Prefilled Syringes and (50) Stat-Flo™ Tips</i>	504640	\$101²⁰

PRODUCT	SKU	25% OFF PRICE
Bulk Syringe Kit <i>Includes: (1) 30mL Prefilled Syringe, (20) Stat-Flo™ Tips, (20) 1.2mL Empty Syringes and (1) Docking Port</i>	504650	\$42¹⁸
Bulk Syringe Refill <i>Includes: (1) 30mL Prefilled Syringe and (1) Docking Port</i>	504660	\$28⁶⁶



Seamfree Composite Wetting Resin & Lubricant

25% OFF

3mL Bottle Kit
Originally \$53⁰⁵
\$39⁷⁹



PRODUCT	SKU	25% OFF PRICE
1 Syringe Kit <i>Includes: (1) 1.2mL Syringe and (10) Exacta-Flo™ Tips</i>	90025	\$15 ⁵³
5 Syringe Kit <i>Includes: (5) 1.2mL Syringes and (50) Exacta-Flo™ Tips</i>	90026	\$74 ⁵⁵
3mL Bottle Kit <i>Includes: (1) 3mL Bottle, (1) Single-Cavity Well and (25) applicators</i>	90027	\$39 ⁷⁹

Designed to make the handling of composite material easier and more consistent. Additionally, **Seamfree Composite Wetting Resin** is unfilled and does not build film thickness or add a visible layer between composite layers.

Seamfree is compatible with all methacrylate materials and can be used in all direct or indirect restorations.

Velvet Etch™

12pk Syringes
Originally \$52⁵³
\$39⁴⁰

Velvet Etch™ is a purple 40% Phosphoric Acid Etch with optimized viscosity for precision placement. Velvet Etch™ is the ideal viscosity, providing accurate placement, control and consistency. It can be applied to vertical tooth surfaces and remain in place and is great for selective etch techniques.



25% OFF

PRODUCT	SKU	25% OFF PRICE
4pk Syringes <i>Includes: (4) 1.2mL Prefilled Syringes and (25) 22ga Pre-Bent Tips</i>	503315	\$17 ⁰⁷
12pk Syringes <i>Includes: (12) 1.2mL Prefilled Syringes and (75) 22ga Pre-Bent Tips</i>	503300	\$39 ⁴⁰