

# Tips on Tips

## Universal Tips

## Tissue Management Tips

## Endodontic Tips

## Periodontal Tips

### Your Total Tip Solutions

Dental syringes are some of the most common pieces of equipment used by dental professionals in a wide range of dental procedures, from the most routine to the most complicated surgical procedures. At Vista Apex, we offer complete syringes, dental syringe tips (including bendable tips, dripless tips, brush tips, and more), and accessories like tip covers.

With many of these products, we offer the option to buy in bulk so your dental practice is always stocked with dental syringes and syringe tips. As with all our dental products, we aim to deliver the highest quality materials to make sure your practice is treating patients right.

SCAN HERE  
For More Information



VISTA APEX  
vistaapex.com  
877.418.4782

**Tissue Management**

**Stat-Flo™ Brush Tip**

- Padded-end brush tip



**20 cup**  
312102

**100 cup**  
312100

**500 bag**  
312105

**Designed for:**  
Hemostatic agents

**MST™ Universal Tip**

- Material saving technology®



**100 cup**  
22ga  
315310

**Designed for:**  
Flowable composites, pulp capping, etchants, & sealants

**Spira-Flo™ Brush Tip**

- Nylon brush fibers for controlled delivery
- Internal spiral
- Chamber minimizes bubbles



**100 cup**  
321400

**500 bag**  
321450

**Designed for:**  
Etchants and resins

**Sol-u-Flo™ Applicator Tip**

- 19ga tip with adjustable brush fiber length
- Curved cannula for optimal entry



**100 cup**  
312200

**500 bag**  
312205

**Exacta-Flo™ Brush Tip**

- Adjustable, soft brush fibers minimize bubbles



**100 cup**  
321610

**500 bag**  
321650

**Designed for:**  
Aqueous solutions such as Caries Indicator, Silane, and Chlorhexidine

**Micro-Flo™ Tip**

- Polyimide tip for microscope enhanced dentistry



**50 cup**  
26ga  
406626

**27.5ga**  
406627

**30ga**  
406630

**Designed for:**  
Flowable composites, pulp capping, etchants, & sealants

**Volu-Flo™ Tip**

- Delivers large volumes of material (18ga)



**100 cup**  
321100

**Designed for:**  
Viscous materials, pastes & putties

**½" Bendable Applicator Tip**

- Blunt end, dead-soft needle resists crimping



**100 cup**  
22ga  
315222

**100 cup**  
25ga  
315225

**Designed for:**  
Etchants, resins, and flowable composites

**Blue-Flo™ Tip**

23ga soft plastic tip




**20 cup**  
316202

**100 bag**  
316222

**Designed for:**  
Irrigation, aeration & impression materials

**Appli-Vac™ – ¾" Bendable Tip**

- Blunt end, dead-soft needle resists crimping



**100 cup**  
16ga  
315116

**100 cup**  
18ga  
315118

**100 cup**  
19ga  
315119

**100 cup**  
20ga  
315120


**100 cup**  
22ga  
315122

**100 cup**  
25ga  
315125

**Designed for:**  
Sealants, flowable composites, etchants, resins, primers

**Pre-Bent Tip**

- Blunt end pre-bent needle
- Opaque photo sensitive hub (22/25ga)



**100 cup**  
16ga  
315316

**100 cup**  
20ga  
315320

**100 cup**  
22ga  
315322

**100 cup**  
25ga  
315325

**Designed for:**  
Etchants, resins and flowable composite

### MST Visco-Tip™

- Bendable, flexible 25ga tip



**20 cup**  
312302  
**50 cup**  
312300

#### Designed for:

Viscous materials such as bioceramic sealers and calcium hydroxide

### Voyager Irrigating Tip

- Dual side ports for safety and efficacy
- Enhanced flexibility to navigate canals
- Silicone stopper for depth reference



|                               |  |                               |                               |
|-------------------------------|--|-------------------------------|-------------------------------|
| <b>17mm</b><br>27ga<br>315727 | <b>17mm Dual Port or Blunt End</b><br>29ga<br>315729 or 315929 | <b>17mm</b><br>30ga<br>315730 | <b>17mm</b><br>31ga<br>315731 |
| <b>27mm</b><br>27ga<br>315827 | <b>27mm</b><br>29ga<br>315829 or 316929                        | <b>27mm</b><br>30ga<br>315830 | <b>27mm</b><br>31ga<br>315831 |

#### Designed for:

Endo irrigating solutions

### Appli-Vac™ 1" Bendable Irrigating Tips

- Slotted and side-vented



|                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|
| <b>20 cup</b><br>23ga<br>315023  | <b>20 cup</b><br>27ga<br>315027  | <b>20 cup</b><br>30ga<br>315030  |
| <b>100 bag</b><br>23ga<br>315523 | <b>100 bag</b><br>27ga<br>315527 | <b>100 bag</b><br>30ga<br>315530 |

#### Designed for:

Endo Irrigating solutions (27 ga / 30 ga)  
Perio Irrigation (23 ga)

### Secure-Lock™ Appli-Vac™ 1" Irrigating Tip

- Slotted and side-vented
- Secure-Lock™ threads lock tip in place




|                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|
| <b>80 bag</b><br>23ga<br>314523 | <b>80 bag</b><br>27ga<br>314527 | <b>80 bag</b><br>30ga<br>314530 |
|---------------------------------|---------------------------------|---------------------------------|

#### Designed for:

Endo Irrigating solutions (27 ga / 30 ga)  
Perio Irrigation (23 ga)

### Vista-Probe™ Irrigating Tip

- Closed-end for side-port delivery
- Lateral delivery of irrigants reduces apical pressure



|                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|
| <b>20 cup</b><br>23ga<br>315123  | <b>20 cup</b><br>27ga<br>315127  | <b>20 cup</b><br>30ga<br>315130  |
| <b>100 bag</b><br>23ga<br>315623 | <b>100 bag</b><br>27ga<br>315627 | <b>100 bag</b><br>30ga<br>315630 |

#### Designed for:

Endo Irrigating solutions (27 ga / 30 ga)  
Perio Irrigation (23 ga)

### Secure Lock™ Vista-Probe™ Tip

- Closed-end for side-port delivery
- Lateral delivery of irrigants reduces apical pressure
- Secure-Lock™ threads lock tip in place



|                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|
| <b>80 bag</b><br>23ga<br>314623 | <b>80 bag</b><br>27ga<br>314627 | <b>80 bag</b><br>30ga<br>314630 |
|---------------------------------|---------------------------------|---------------------------------|

#### Designed for:

Endo Irrigating solutions (27 ga / 30 ga)  
Perio Irrigation (23 ga)

### Elasti-Vac™ Tip

- Narrow, flexible tip easily navigates curved canals



**50 cup**  
406150

#### Designed for:

Endodontic evacuation

### Versa-Brush™

- Rotational utility brush



**20 cup**  
204320  
**40 cup**  
204300

#### Designed for:

Cleaning post holes and coronal access openings

### NiTi SuperFlex™ Tip

- AUTOCLAVABLE. Slotted and side-vented, 30ga for safe irrigation
- 2.5x more flexible



|   |  |                         |
|---|--|-------------------------|
| <b>17mm (short)</b><br><b>6 cup</b><br>405107 | <b>25mm (long)</b><br><b>6 cup</b><br>405108 | <b>20 cup</b><br>406320 |
| <b>Mixed 6 cup</b><br>405106                  | <b>50 cup</b><br>406350                      |                         |

#### Designed for:

Endo irrigating solutions

### Flexi-Glide™ Utility Tip

- 30ga flexible tip for easy access



**20 cup**  
406320  
**50 cup**  
406350

### Periodontal

#### Pocket-Probe™

- Closed-end with multiple perforations for thorough flushing



**100 cup**  
21ga  
316102

#### Designed for:

Flushing of perio pockets or inter-dental spaces