




Material Saving Tip

Sometimes it's the Little Things that Make all of the Difference

- The **MST Tip** minimizes the internal volume of the tip greatly **reducing waste, up to 87%**!
- Ideal solution to get the best value from costly materials such as bioceramic sealers.
- Unique, flexible tip improves flow and placement of viscous materials.
- Rotatable and bendable shoulder for optimum positioning with non-clogging design.




ENDODONTIC TIPS

MST Visco-Tips™

- Bendable, flexible 25ga tip
- Easily trimmable to shorter lengths
- Improves flow and placement of viscous materials

Designed for: Bioceramic sealers, and Calcium Hydroxide

20 pk 312302
50 pk 312300



UNIVERSAL TIPS

MST Tips™

- 22ga pre-bent tip
- Opaque photo sensitive hub

Designed for: Costly flowable composites, pulp capping materials, etchants, sealants

100 pk 315310



NEW!

ENDODONTIC TIPS


MST Visco-Tips™ Stainless Steel

- Annealed, 25ga stainless steel tip that is easily bendable
- Provides better tactile feedback during canal navigation
- Available in 17mm and 22mm working lengths

Designed for: Bioceramic sealers, and Calcium Hydroxide

50 pks 17mm 312411 22mm 312412

Use more of what you paid for!




Conventional tips trap a significant amount of unusable material inside of the hub. The structure of the MST Tip and the bendable, flexibility allows majority of the material to go through.

We Did the Math on Savings

	Leading Competitor	MST Visco Tip
Procedures / syringe	8	37
Endodontists procedures / day	5.7	
Endodontists procedures / year (200 working days)	1140	
Number of syringes needed / year	143	31
Cost per syringe	\$185	
Annual sealer cost	\$26,455	\$5,735
SAVINGS USING MST	\$20,720	78%

Save Material Save Money



Save money and material on tips with Vista Apex! From applicator tips to air/water tips, Vista Apex has it all.

VA VISTA APEX
ENDO RESTO LET'S GO!

