


# ONE TIP FOR ALL SOLUTIONS

Enhance Your Endo  
Irrigation with the V.I.P. Tip



**Our New V.I.P. Tip  
(Voyager Irrigating  
Precision Tip) is one  
tip for all solutions.  
This Universal Irrigating  
Tip can be used for  
all Endo Irrigating  
Procedures**

Pre-bent at the optimal  
angle for access, easily  
adjusted by hand or  
with pliers.



Small flexible cannula  
to navigate canals  
and deliver irrigants  
or sealers.

## What makes a V.I.P. Tip Better?


	Gauges Are Dual Side Port	Available in Multiple Lengths	Enhanced Flexibility	Low Cost
V.I.P. Tip	✓	✓	✓	✓
Ultradent® NaviTip™	✗	✓	✓	✗
PacDent OptiProbe™	✗	✓	✗	✓

*Ultradent® and NaviTip™ are registered trademarks of ULTRADENT® Products, Inc.  
OptiProbe™ is a registered trademark of PacDent, Inc.*

## Spotlight Performance

One tip for all solutions! Enhance your Endo Irrigation with the V.I.P. Tip (Voyager Irrigating Precision Tip).

To ensure safe delivery, this tip was engineered so that irrigants are expressed out the sides of the tip. Using a metal shoulder in place of a large plastic hub, the two lengths of the V.I.P. Tip offers the same full range of working lengths compared to competitors without having to stock as many variations in your operator. The reduced bulk of the hub also allows for easier access in smaller oral openings.



## Available in Two Lengths, 17mm and 27mm

GA	STYLE	LENGTH	COLOR	PART #
29GA	Dual Side Port	17mm 27mm	Orange	315729 315829

- Dual side ports for safety and efficacy
- Enhanced flexibility to navigate canals
- Available in two lengths to help fit in limited openings reach all canal depths
- Silicone stopper for depth reference
- Pre-bent at the optimal angle for access, easily adjusted by hand or with pliers
- Stainless steel shoulder for zero-flow impedance and maximum visualization
- Secure-Lock hub for reliable syringe attachment

## Order Now!

Call toll free 1-877-418-4782 or visit [vistaapex.com](http://vistaapex.com)