

# Seal. Repair. Regenerate.



## RE-GEN™ Bioactive EndoPutty

RE-GEN™ Bioactive EndoPutty is a bioactive bioceramic root and pulp treatment designed for permanent root repair, including perforation, resorption, apexification and pulp capping.

Its advanced calcium sodium phosphosilicate and calcium aluminate formulation ensures superior washout resistance, reliable hermetic seal, and excellent handling characteristics.



## RE-GEN™ Bioactive EndoSealer

RE-GEN™ Bioactive EndoSealer is a hydraulic, ready-to-use injectable bioactive paste that utilizes moisture naturally present in dentinal tubules to initiate a setting reaction. Its excellent physical and chemical properties make it the obvious root canal sealer.



*"I absolutely love RE-GEN™ Bioactive EndoPutty – great price, quick setting time and a packaging design that prevents the material from setting in the container. RE-GEN™ Bioactive EndoPutty is my go-to choice for perforation repairs, pulpotomies, apical surgeries and all other endodontic procedures."*

— Dr. Randall Cross

Fast  
& Ready  
to Use

## A powerful combination

Achieve a fluid tight seal throughout the canal with RE-GEN™ Bioactive EndoSealer prior to repairing with RE-GEN™ Bioactive EndoPutty. These bioactive combination products promote bio-regenerative healing due to their Bioglass and BioCeramic components.

Choosing a combination that is effective, fast and ready to use is not only a good decision – it's an economical one, too.

### RE-GEN™ Bioactive EndoSealer



#### Key benefits:

- **Ready to use** – Single component, highly radiopaque bioactive paste
- **High pH** – Creates high alkaline environment to prevent bacteria proliferation
- **Fast-setting** – Efficient set time of under 1 hour
- **Dimensionally stable** – Low solubility with high wash out resistance
- **Cost-effective** – Paired with material saving technology MST Visco-Tip™ for low procedure cost

#### Bioactive EndoSealer Kit Options\*

504753	<b>FULL KIT</b> <ul style="list-style-type: none"> <li>• 1 x 2g prefilled syringe</li> <li>• 20 x MST Visco-Tips</li> </ul>
50474	<b>FULL KIT</b> <ul style="list-style-type: none"> <li>• 1 x 2g prefilled syringe</li> <li>• 20 x MST Visco-Tips stainless steel</li> </ul>

### RE-GEN™ Bioactive EndoPutty



#### Key benefits:

- **Ready to use** – A single-component, easy-to-use putty that's highly radiopaque
- **Stays in place** - Resists washout after
- **Fast-setting** – Cures in under an hour, enabling quicker treatments
- **Antibacterial properties** – High pH helps inhibit bacterial growth
- **Modable and versatile** - Firm yet adaptable to fit various clinical needs and anatomies

#### Bioactive EndoPutty\*

504756	<ul style="list-style-type: none"> <li>• 1 x 1g prefilled syringe</li> </ul>
--------	--

\* Product has a shelf life of 24 months from manufacturing date.

#### Vista Apex Solutions

2200 South Street  
Racine, WI 53404

Phone: (877) 857-8213

Email: [info@vista-dental.com](mailto:info@vista-dental.com)

Website: <https://vistaapex.com>