

Tips on Tips

Your Total Tip Solutions

VistaApex Solutions provides complete dental syringes, a wide range of syringe tips (bendable, dripless, brush, and more), and accessories like tip covers. Many items are available in bulk to keep your practice stocked. Our products are made with the highest quality materials to support excellent patient care.



SCAN HERE
to order now or for
more information

PRE-BENT DELIVERY

Pre-Bent Tip

· Opaque photo sensitive hub (20/22ga)

Micro-Flo™ Tip

· Polyimide tip for microscopic enhanced dentistry



100 cup
315316



100 cup
315320



100 cup
MST
315310
Standard
315322



100 cup
315325



50 cup
26ga 406626
27.5ga 406627

16ga

20ga

22ga

25ga

26ga & 27.5ga

Designed for:

Etchants, resins, flowable composite, resin barriers, in-office whitening

SPECIALTY DELIVERY TIPS

Spira-Flo™ Brush Tip

· Nylon brush fibers control delivery
· Internal spiral minimizes bubbles



100 cup
321400
500 bag
321450

18ga

Designed for:

Etchants, Bonding Agents, Flowable Composites, Pit & Fissure Sealants

Volu-Flo™ Tip

· Delivers large volumes of material



100 cup
321100

18ga

Designed for:

Viscous materials, pastes & putties

Exacta-Flo™ Brush Tip

· Adjustable, soft brush fibers minimize bubbles



100 cup
321610
500 bag
321650

18ga

Designed for:

Aqueous solutions such as Caries Indicator, Silane, Chlorhexidine, Primers and Resin Coats

Sol-u-Flo™ Applicator Tip

· Adjustable brush length
· Curved cannula for optimal entry



100 cup
312200
500 bag
312205

19ga

SPECIALTY DELIVERY TIPS

Stat-Flo™ Brush Tip

· Padded-end brush tip



20 cup
312102
100 cup
312100
500 bag
312105

19ga

Designed for:

Hemostatic agents, bonding agents, iron solutions

Blue-Flo™ Tip

· Soft plastic tip



20 cup
316202
100 bag
316222

23ga

Designed for:

Irrigation, aeration, impression materials, caulks, putties, scrubs and slurries

½" Bendable Applicator Tip

· Dead-soft needle resists crimping



100 cup
22ga 315222
25ga 315225

22ga & 25ga

Designed for:

Etchants, Bonding Agents, Flowable Composites, Pit & Fissure Sealants

BENDABLE DELIVERY

Appli-Vac™ – ¾" Bendable Tip

· Dead-soft needle resists crimping



100 cup
315116



100 cup
315118



100 cup
315119



100 cup
315120



100 cup
315122



100 cup
315125

16ga

18ga

19ga

20ga

22ga

25ga

Designed for:

Sealants, flowable composite, etchants, resins, primers. Ideal tip for application and aspiration. Use with Luer HVE Untipped Adapter for removal of moisture from endodontic canals (315250)

IRRIGATION

Appli-Vac™ 1" Bendable Irrigating Tip

- Slotted and side-vented for safe irrigation



23ga Standard
20 cup 315023
100 cup 315523



27ga Standard
20 cup 315027
100 cup 315527



30ga Standard
20 cup 315030
100 cup 315530

23ga Secure Lock
80 bag 314523

27ga Secure Lock
80 bag 314527

30ga Secure Lock
80 bag 314530

23ga

27ga

30ga



20 cup
23ga 315123
27ga 315127
30ga 315130



100 bag
23ga 315623
27ga 315627
30ga 315630

23ga, 27ga, & 30ga

Designed for:

Endo Irrigating and citric acid solutions (27ga & 30ga) Perio Irrigation (23ga)

NiTi SuperFlex™ Irrigating Tips Autoclavable

- Slotted and side-vented



6 cup
17mm Short
405107
25mm Long
405108

30ga

Designed for:

Endo Irrigating Solutions

IRRIGATION

FEATURED TIP

VIP Tip (Voyager Irrigating Precision Tip)

One tip for all solutions!

- Dual Side Port
- 2 length options



50 cup
17mm 315729



50 cup
27mm 315829



Voyager Irrigating Tip

- Dual Side Port
- 3 gauge and 2 length options



50 cup
17mm
27ga 315727
30ga 315730
31ga 315731



50 cup
27mm
27ga 315827
30ga 315830
31ga 315831

29ga

29ga

27ga, 30ga, & 31ga



SCAN HERE

To download
the Digital
Endo Protocol

Designed for:

Endo Irrigating and citric acid solutions
Use with Endo Sealer Solvent™

EVACUATION

Elasti-Vac™

- Narrow/flexible tip easily navigates curved canals



Blue Straight Tip
50 cup
406150



Clear Pre-Bent Tip
50 cup
406151

25ga

25ga

Designed for:

Endodontic evacuation, irrigation or peri-implant pockets with CHX gel. Use with Luer HVE Untipped Adapter for removal of moisture from endodontic canals (315250)
Tip Diameter: 0.25mm ID, 0.52mm OD

DELIVERY

MST Visco-Tip™ Polymer Trimmable



Polymer
20 cup
312302
50 cup
312300

MST Visco-Tip™ Stainless Steel Bendable



Stainless Steel
17mm
50 cup
312411



Stainless Steel
22mm
50 cup
312412

25ga

25ga

25ga

Designed for:

Viscous materials such as bioceramic sealers, calcium hydroxide, EDTA file lubricant, MTA Cement

Voyager Tip

- Blunt end
- 2 length options



50 cup
17mm 315929
27mm 316929

Blunt 29ga

Designed for:

Delivery of high-flow sealer or paste materials. May also be used for evacuation.

PERIODONTAL

Pocket-Probe™

- Closed-end with multiple perforations for thorough flushing



100 cup
316102

21ga

Designed for:

Flushing of perio pockets or inter-dental spaces. Ideal tool for at home perio irrigation.