

One Tip for all Solutions

Enhance Your Endo Irrigation with the VIP Tip

The VIP Tip (Voyager Irrigating Precision Tip) simplifies endodontic irrigation by using side-expressing ports and a slim metal shoulder for easier access. Available in two lengths, it eliminates the need for multiple tip variations and improves usability in tight spaces.



VIP Tip

- **Easy Access to Canals**
Dual side ports and two lengths make it ideal for safe and effective irrigation—even in hard-to-reach canals.
- **Precise Every Time**
Pre-bent at the optimal angle for access with adjustable silicone stopper helps you reach the right depth with confidence and control.
- **Smooth, Steady Flow**
The Secure-Lock hub and slim stainless-steel shoulder provide a tight syringe fit and clear visibility—no flow blockages.



VIP Tip Options

315729	29GA 17MM • 1 x 50 tips cup
315829	29GA 27MM • 1 x 50 tips cup

VistaApex Solutions

2200 South Street
Racine, WI 53404

Phone: (877) 857-8213
Email: info@vista-dental.com
Website: <https://vistaapex.com>



V+ VistaApex
SOLUTIONS
Value. Added.