

More than a cement, a regenerative solution.

RE-GEN® BIOACTIVE Cement




RE-GEN® Bioactive Cement is a dual-cure, self-adhesive resin cement with built-in bioactivity for indirect restorations. It delivers strong, reliable bonding while actively supporting remineralization at restoration margins and preventing recurrent decay.

Powered by Bioactive technology

RE-GEN® Bioactive Cement goes beyond traditional cementation by incorporating Bioglass technology that actively interacts with the oral environment. Through ion exchange, it supports remineralization and helps create a protective environment at restoration margins—where recurrent decay and breakdown can begin. This bioactive response works alongside strong, reliable adhesion to help support long-term restoration integrity and clinical confidence.

RE-GEN® Bioactive Cement uses a streamlined shade system with three options: **A2**, **White Opaque**, and **Translucent**—to accommodate any clinical situation.



RE-GEN® Bioactive Cement

- | | |
|-------|--|
| 90920 | Shade A2 <ul style="list-style-type: none"> • 1 x 10g Dual-Barrel A2 Syringe • 8 x Endo Intraoral Tips • 8 x Brown Base/Clear Mixing Tip • 1 x Sealable Foil Pouch • 1 x IFU |
| 90921 | Shade White Opaque <ul style="list-style-type: none"> • 1 x 10g Dual-Barrel WO Syringe • 8 x Endo Intraoral Tips • 8 x Brown Base/Clear Mixing Tip • 1 x Sealable Foil Pouch • 1 x IFU |
| 90922 | Shade Translucent <ul style="list-style-type: none"> • 1 x 10g Dual-Barrel TR Syringe • 8 x Endo Intraoral Tips • 8 x Brown Base/Clear Mixing Tip • 1 x Sealable Foil Pouch • 1 x IFU |

RE-GEN® Bioactive Cement Key Features:

- **Bioactive Protection with Bioglass**
Supports remineralization and helps protect margins for long-term restoration integrity, and supports the prevention of recurrent decay.
- **Strong, Reliable Bond Strength**
Delivers durable adhesion across multiple substrates for long-term clinical confidence
- **Fast Self-Cure Performance**
Efficient 3-5 minute self-cure reduces chair time and improves productivity
- **Self-Adhesive Simplicity**
Eliminates the need for separate bonding steps, streamlining your workflow



VistaApex Solutions

2200 South Street
Racine, WI 53404

Phone: (877) 857-8213
Email: info@vista-dental.com
Website: <https://vistaapex.com>

Scan to
learn more

